




di Ing. Sergio Serrano & C.

### ELECTROMEDICAL EQUIPMENT – NATURAL PRODUCTS

VIA RONCHI, 30 – 20134 – MILAN (I) - tel. +39 02 21597558 - fax +39 02 21597559 - E. MAIL: zener@zener-seab.it



### BFBC mod. ZTDS

BFBC is a computerized diagnostic instrument which tracks measurements, according to the techniques of potential recalls (bio feed-back), in peripheral electrodermics.

The interactive software applies EVA, RYODORAKU, VEGA techniques and results in a highly accurate valuation of the oxidative stage in the Organism (homotoxicological stage)

BFBC is a truly essential device when it comes to using PHYTOBIOENERGETIC ZENER to drain and purify the Organism and is propaedeutic in Academic and Natural Medical treatment (Homeopathy, Homotoxocology, Floritherapy, Acupuncture etc.)

BFBC can be applied to carry out biophotonic tests of sample vials simulating the Organism's reaction to the real intake of homeopathic, pharmaceutical and food remedies.

BFBC plays a role of primary importance in studies carried out on food intolerance, allergies and dental equipment and comes with a set of test vials (or filters) which help to make the diagnostic valuation.

BFBC picks up signals resulting from digital testing using oscillators or DCs.

BFBC can be connected to ISOQUANT III to compile ECCIMERI® of ground or digital signals tested by BFBC itself..

#### TECHNICAL FEATURES BFBC

- analogical setting device
- setting knob
- incorporate capacitive winze for vial tests
- push rod with pressure and humidity dial
- earth electrode
- rechargeable battery 24 Volt
- battery recharger included
- 9-pin plug to connect serial computers
- digital signal test input
- output for ISOQUANT III
- size: 320x210x75 mm
- weight : 2,3 Kg

#### TECHNICAL FEATURES COMPUTER MINIMUM

- Operating system Windows 98 or superior
- Microprocessor Pentium 500 Mhz.